

CULTURE 

Roots

TISSUE CULTURE MEDIA
Rooting & Production

 ATHENA®

Sachet makes
125 mL

CULTURE 

Shoots

TISSUE CULTURE MEDIA
Storage & Multiplication


 ATHENA®

Sachet makes
125 mL

ALL INCLUSIVE MEDIA
FORMULA. JUST ADD WATER.

PREPARATION STEPS:

- 1** Pour entire sachet contents into 250 mL mixing vessel.
- 2** Fill vessel to 125 mL with clean water.
- 3** Close and agitate until dissolved.

 **WARNING: KEEP OUT OF REACH OF CHILDREN & ANIMALS. NOT FOR HUMAN CONSUMPTION. FOR LAB USE ONLY.**